

D.I.S

D.I.S Germany GmbH

BREITE STRASSE 2

01796 PIRNA

GERMANY

Tutorial IBT-001

Introduction of Materials into Ion Sources

Use of Volatile Chemical Compounds

Author[©] : Günter Herbert Zschornack

December 6, 2021

The use of metal ions is required for a wide range of applications. Metal ions are often difficult to introduce into ion sources and known techniques such as evaporation furnaces or sputtering techniques do not always lead to the desired success. Therefore, the MIVOC¹ method for generating metal ions has proven to be favorable for a large number of metals and was first proposed by [Ärje et al. 1993]. This method uses metallo-organic compounds whose vapor pressures at room temperature allow the compound to be introduced as a gas into the ion source via a precision inlet valve.




Figure 1: Principle diagram for the MIVOC method for loading the ionization volume of ion sources. p1 – saturated vapour pressure of the used chemical compound, p2 – pressure inside the ion source, c – conductance of the transmission line, v1 – valve for a pump connection, v2 – lock valve.

The generalized principle scheme for the use of the MIVOC method is shown in Figure 1 for a typical application with an ECRIS. The substance used is placed in a MIVOC chamber and a saturated vapor pressure is formed at a certain vacuum with the pressure p1. To ensure clean conditions when introducing the gas into the ion source, the chamber area and the gas transport section are first pumped down to the pressure p1. Connection to a backing pump is made via valve v1 with valve v2 closed. After the appropriate working pressure has been reached, the pumping process is terminated by closing the valve v1 in the direction of the pump, and valve v2 is opened. The gas can then flow into the ion source at pressure p2. The design of the gas transport section with a conductance c based on this must also be taken into account here.

Since even small amounts of gas are sufficient for source loading, it is possible to operate with compounds of metals, solids or liquids that have a vapour pressure of the order 10^{-3} mbar (for example for ECRIS) or lower of up to about 10^{-5} mbar for EBIS. The limit for the lowest usable vapour pressure here is the conductance of the gas injection system.

¹MIVOC: Metal Ions from Volatile Compounds

The technical realization of the MIVOC method can generally be done in two ways:

1. If the vapor pressure at room temperature is about 10 mbar or more, a mass flow controller with a piezoelectric valve can be used.
2. Vapor pressures of about 0.01 mbar are realized via the temperature control of the MIVOC chamber and the material inlet into the ion source.

To avoid condensation in the gas transfer line, it is necessary to control the temperature profile in the pipeline. A temperature higher than that of the MIVOC chamber is maintained in the gas transfer system. This can be achieved by means of corresponding, possibly regulated heating coils in the gas transfer system.

Values for vapour pressures of various compounds are often not available in the literature. This makes it necessary to calculate vapor pressures from available data sets. If the temperature of the melting point or the entropy of vaporization at the melting point are known, the vapour pressure at room temperature can be estimated [Koivisto et al. 1994, Ullrich 1997, Ullrich et al. 1998].

Starting from the Clausius-Clapeyron equation which is approximately valid for vapor pressure curves

$$\frac{d \ln P_{V_p}}{dT} = \frac{\Delta S_V}{\Delta Z RT^2} \quad (1)$$

with P_{V_p} – vapor pressure, T – temperature, ΔS_V – enthalpy of vaporization, ΔZ – compressibility and R – molar gas constant the following simplified relationship is obtained [Weast 1978]

$$\ln P_{V_p} [\text{mbar}] \approx 0.75 \frac{\Delta S_v [\text{kJ/mol}] \Delta T [\text{K}]}{0.15 \Delta T [\text{K}] - T [\text{K}]} + 2.881 \quad (2)$$

ΔS_v is the entropy of vaporization² and is generally available in tabular form. Further applies

$$\Delta T = T - T_b$$

with T_b – boiling point of the compound.

For a more precise estimation, a modified Watson method [Lyman et al. 1990] can be used, which can be derived from the Clausius-Clapeyron equation (1) and whose range of application includes pressures between 10^{-3} mbar to about 10^{-7} mbar. The method described here only requires knowledge of the melting temperature T_b as substance-specific parameter. A more detailed description of the modified Watson method is given by [Ullrich 1997] which is followed here.

²The vaporization entropy ΔS_v describes the change in entropy of a liquid during the transition to the gas phase.

The temperature dependence of the evaporation enthalpy can be determined via

$$\Delta S_V = \Delta S_{V_b} \left(\frac{1 - T/T_c}{1 - T_b/T_c} \right)^m \quad (3)$$

with ΔS_{V_b} – enthalpy of vaporization at the melting point, T_c – critical temperature³ and m – substance-specific parameter.

Here for liquids applies $m = 0.19$ and for solids yields

$$m = \begin{cases} 0.36 & \text{for } T_{Q_b} > 0.6 \\ 0.80 & \text{for } 0.6 > T_{Q_b} > 0.5 \\ 1.19 & \text{for } T_{Q_b} < 0.5 \end{cases}$$

with $T_{Q_b} = T/T_b$. Assuming

$$T_C \approx \frac{3}{2} T_b$$

Equation (3) can be transformed to

$$\Delta S_V \approx \Delta S_{V_b} (3 - 2T_{Q_b})^m \quad (4)$$

Inserting Equation (4) in equation (1) gives

$$\ln P_{V_p} \approx - \frac{\Delta S_{V_b}}{\Delta Z_b R T_b} \left[\left| \frac{(3 - 2T_{Q_b})^m}{T_{Q_b}} + 2m (3 - 2T_{Q_b})^{m-1} \ln T_{Q_b} \right|_1^{T_{Q_b}} \right. \\ \left. - 4m(m-1) \underbrace{\int_1^{T_{Q_b}} (3 - 2T_{Q_b})^{m-2} \ln T_{Q_b} dT_{Q_b}}_{\approx 0} \right]$$

and thus final

$$\ln P_{V_p} \approx \frac{\Delta S_{V_b}}{\Delta Z_b R T_b} \left[1 - \frac{(3 - 2T_{Q_b})^m}{T_{Q_b}} + 2m(3 - 2T_{Q_b})^{m-1} \ln T_{Q_b} \right] \quad (5)$$

³The critical temperature is the temperature of substances at which the liquid-gas phase transition proceeds without any change in volume. As the temperature belonging to the critical point, it defines the end point of the vapor pressure curve in the phase diagram and is thus the temperature above which liquefaction of a gas is no longer possible, no matter how high the pressure.

According to [Lyman et al. 1990] the entropy at the melting point can be described via

$$\frac{\Delta S_{V_b}}{T_b} = K_S (8.75 + R \ln T_b) \quad (6)$$

with $R = 1.987 \text{ cal}/(\text{mol}\cdot\text{K})$ and T_b in Kelvin. K_S is a tabulated structure coefficient.

If equation (6) is substituted into equation (5), the vapor pressure can be calculated in atmospheres or in multiples of 760 Torr. The compressibility can be approximated to $\Delta Z_b = 0.97$ in this estimation [Lyman et al. 1990]. [Ullrich 1997] indicates that the error of the calculations is below 50% in the pressure range between 10^{-3} mbar up to 10^{-7} mbar.

It should be noted that the saturation evaporation pressure is constant at constant temperature. The time it takes for the saturation evaporation pressure to build up in a closed volume can be calculated according to [Koivisto et al. 1994]

$$t_{sat} [\text{s}] = \frac{p_{sat} [\text{mbar}] V_{MIVOC} [\text{l}]}{N_{max} \left[\frac{\text{mbar l}}{\text{s mg}} \right] m [\text{mg}]} \quad (7)$$

with p_{sat} – saturation pressure, V_{MIVOC} – volume of the MIVOC chamber, N_{max} – maximum evaporation rate of the compound and m – mass of the used compound. With commonly used MIVOC chambers, the saturation time is a maximum of a few 10 ms. It is also important to know that with common MIVOC chambers, their volume does not affect the partial pressure of the compound used. Small chambers can therefore be used, where moisture can be more easily removed from the chamber.

Table 1 compiles suitable metallo-organic compounds for use in MIVOC technology. When using the MIVOC method, it should be noted that a number of the reagents may be toxic, corrosive, hygroscopic or have other important properties. Therefore, it is recommended to always inform yourself about the recommended handling precautions of various substances before using them.

It should also be mentioned that the breakup of the molecules used produces larger quantities of hydrogen and carbon ions, which makes it difficult to set specific ion charge distributions in the ion source. Carbon also contaminates the source and thus leads to a reduction in the ion currents of the target element.

Table 1: Organometallic and other compounds for loading ion sources with metal ions. The table is based on a compilation by [Ullmann 2005], including various additions by the author. The specification of a compound without a literature reference means that the listed compound can potentially be used, but that the author is not aware of any literature work on it.

Z	name	sum formula	references
5	1,2-dicarbado-decaborane(12) also called decaborane(12)	$C_2H_{12}B_{10}$	[Bogomolov et al. 1999] [Koivisto et al. 2001] [Bogomolov et al. 2015] [Loginov et al. 2018]
	decaborane(14)	$B_{10}H_{14}$	[Oshiro et al. 2014]
12	magnesocene	$Mg(C_5H_5)_2$	[Koivisto et al. 1998] [Takasugi et al. 2010]
13	trimethyl aluminium	AlC_3H_9	
14	silicon tetrachloride	$SiCl_4$	[Takasugi et al. 2010]
	tetramethylsilane	$Si(CH_3)_4$	[Takasugi et al. 2010] [Bondarchenko et al. 2020]
	tetrakis(trimethyl-silyl)silane	$Si[(CH_3)_3Si]_4$	[Koivisto et al. 2002]
15	diethyl phosphite	$(C_2H_5O)_2P(O)H$	[Kentsch et al. 2004a]
22	titanium (IV) chloride	$TiCl_4$	[Kentsch et al. 2004] [Takasugi et al. 2010]
	titanium tetrafluoride	TiF_4	[Bogomolov et al 2014]
	cyclopentadienyl cycloheptatrienyl titanium	$C_5H_5TiC_7H_7$	[Bogomolov et al 2014]
	(trimethyl)pentamethyl cyclopentadienyl titanium(IV)	$(CH_3)_5C_5Ti(CH_3)_3$	[Bogomolov et al 2014] [Bondarchenko et al. 2020] [Pugachev et al. 2021]
	titanocene dichloride	$(C_5H_5)_2TiCl_2$	
	titanium (IV) ethylate	$Ti(OC_2H_5)_4$	
23	vanadocene	$V(C_5H_5)_2$	[Ullrich et al. 1998] [Bogomolov et al. 2015] [Bondarchenko et al. 2016]
24	chromozene	$Cr(C_5H_5)_2$	[Koivisto et al. 1998] [Bogomolov et al. 2015] [Bondarchenko et al. 2016] [Bondarchenko et al. 2020] [Pugachev et al. 2021]

Z	name	sum formula	references
	chromium hexacarbonyl	$\text{Cr}(\text{CO})_6$	[Koivisto et al. 1998] [Nakagawa et al. 1998]
25	methylcyclopentadienyl mangan-tricarbonyl	$\text{Mn}(\text{CH})_5\text{CH}_3(\text{CO})_3$	[Werner et al. 2000] [Kentsch et al. 2002a]
26	ferrocene	$\text{Fe}(\text{C}_5\text{H}_5)_2$	[Koivisto et al. 1994] [Stiebing et al. 1999] [Werner et al. 2000] [Kentsch et al. 2002] [Ullmann 2005] [Takasugi et al. 2010] [Bondarchenko et al. 2020] [Liu et al 2020]
27	cobaltocene	$\text{Co}(\text{C}_5\text{H}_5)_2$	[Koivisto et al. 1998] [Bogomolov et al. 2015] [Bondarchenko et al. 2016] [Bondarchenko et al. 2020] [Loginov et al. 2020]
	dicobalt octacarbonyl	$\text{Co}_2(\text{CO})_8$	[Koivisto et al. 2002]
	cyclopentadienyl-cobalt- dicarbonyl	$\text{CoC}_5\text{H}_5(\text{CO})_2$	
28	nickelocene	$\text{Ni}(\text{C}_5\text{H}_5)_2$	[Koivisto et al. 1994] [Stiebing et al. 1999] [Barue et al. 2000] [Kentsch et al. 2002a] [Kentsch et al. 2002] [Bogomolov et al. 2015] [Bondarchenko et al. 2016] [Bondarchenko et al. 2020] [Loginov et al. 2020]
29	copper (II) acetyl acetate	$\text{Cu}(\text{CH}_3\text{COCHCOCH}_3)_2$	
30	diethylzinc	$\text{Zn}(\text{C}_2\text{H}_5)_2$	[Kentsch et al. 2004]
	4-cyclohexylbutyric acid zinc salt	$\text{Zn}[\text{C}_6\text{H}_{11}(\text{CH}_2)_3\text{CO}_2]_2$	
31	trimethylgallium	GaC_3H_9	
	gallium (III) chloride	GaCl_3	

Z	name	sum formula	references
32	tetramethylgermanium	$\text{Ge}(\text{CH}_3)_4$	[Werner et al. 2000] [Bogomolov et al. 2015] [Bondarchenko et al. 2016] [Bondarchenko et al. 2020]
	tetramethylgermane	$\text{Ge}(\text{CH}_2\text{CH}_3)_4$	[Bogomolov et al. 2015] [Bondarchenko et al. 2016]
	dimethylgermane	$\text{Ge}(\text{CH}_3)_2\text{H}_2$	[Takasugi et al. 2010]
33	phenylarsine	$\text{As}(\text{C}_6\text{H}_5)\text{H}_2$	
	arsenic trioxide ethyl ester	$\text{As}(\text{C}_2\text{H}_5\text{O})_3$	
34	dimethylselenide	SeC_2H_6	
39	yttrium (III) isopropylate	$\text{Y}[(\text{CH}_3)_2\text{CHO}]_3$	
40	zirconocene dichloride	$(\text{C}_5\text{H}_5)_2\text{ZrCl}_2$	
41	niobium (V) ethylate	$\text{Nb}(\text{OC}_2\text{H}_5)_5$	
42	molybdenum hexacarbonyl	$\text{Mo}(\text{CO})_6$	[Koivisto et al. 1998] [Nakagawa et al. 1998]
44	ruthenocene	$\text{Ru}(\text{C}_5\text{H}_5)_2$	[Koivisto et al. 1998] [Nakagawa et al. 1998]
48	dimethylcadmium	$\text{C}_2\text{H}_6\text{Cd}$	
49	triethylindium	$(\text{C}_2\text{H}_5)_3\text{In}$	
	indium (III) tert-butylate	$[(\text{CH}_3)_3\text{CO}]_3\text{In}$	
50	Tetramethyltin	$(\text{CH}_3)_4\text{Sn}$	[Werner et al. 2000]
51	Trimethylantimon	$\text{C}_3\text{H}_9\text{Sb}$	[Nakamura et al. 2000]
	antimony (III) chloride	SbCl_3	
	triisopropyl antimony acid ester	$\text{Sb}[\text{OCH}(\text{CH}_3)_2]_3$	
52	dimethyltellurium	$\text{C}_2\text{H}_6\text{Te}$	[Nakamura et al. 2000]
53	diiodomethane	I_2CH_2	[Koivisto et al. 2002]
	diiodine	I_2	
59	tris(cyclopentadienyl) praseodymium	$(\text{C}_5\text{H}_5)_3\text{Pr}$	
60	neodymium (III) cyclo-pentadienylide	$(\text{C}_5\text{H}_5)_3\text{Nd}$	
65	terbium (III) acetyl acetate hydrate	$[\text{CH}_3\text{COCH}=\text{C}(\text{O}-)\text{C}(\text{CH}_3)_3]_3\text{Tb}\cdot x\text{H}_2\text{O}$	

Z	name	sum formula	references
72	bis-(tert.-butylcyclopentadienyl)dimethyl hafnium	$C_{18}H_{26}Cl_2Hf$	
	dimethylbis(cyclopentadienyl)hafnium(IV)	$(C_5H_5)_2Hf(CH_3)_2$	[Jovovic et al. 2007] [Bogomolov et al. 2015] [Bondarchenko et al. 2016] [Bondarchenko et al. 2020]
74	tungsten hexacarbonyl	$W(CO)_6$	[ZiMin et al. 2000] [Koivisto et al. 2002]
76	osmocene	$Os(C_5H_5)_2$	[Koivisto et al. 1998] [Nakagawa et al. 1998] [Koivisto et al. 2002]
78	(cis,cis-1,5-Cyclooctadien)-dimethylplatin (II)	$C_{10}H_{18}Pt$	
83	bismuth (III) chloride	$BiCl_3$	
	triphenyl bismuth	$(C_6H_5)_3Bi$	
	bismuth(III)-2,2,6,6-tetramethyl-3,5-heptanedionate	$[(CH_3)_3CCOCH=C(O-)C(CH_3)_3]_3Bi$	

In combination with the MIVOC method, ion currents for charge state separated metal ions up to approx. some $100\ \mu A$ were achieved in ECRIS. For EBIS, DC currents up to a maximum of nA can be achieved for low ion charge states, for higher charged ions significantly less (about pA).

As examples, ion beam currents are tabulated in Table 2 as they were published in connection with the operation of EBIS and ECRIS. In the case of ECRIS, ion currents are specified in the DC mode, whereas pulsed or DC currents are specified for EBIS. The values measured for EBIS as DC currents were obtained in leaky mode, i.e. from a continuous ion current across the electrostatic trap potential.

Table 2: Ion currents during loading of ECRIS and EBIS with organometallic and other compounds.

For the sake of compactness, the individual ECRIS and EBIS used are not characterized in detail and the operating parameters sets are not listed. However, the interested reader will find these in the original literature cited for each value. The specified ion currents are to be understood as electric particle currents, i.e. to determine the extracted ion numbers, the specified values must still be divided by the indicated ion charge state. Furthermore, it should be noted that the tabulated values are only orientation values. In some cases, the values were obtained from spectra in which the maximum extractable current was optimized for a specific ion charge state. Careful optimization can therefore lead to even higher ion currents for individual ion charge states. In the case of EBIS, "*" is used to characterize ion currents measured in a pulsed EBIS regime. The abbreviation "Im" stands for leaky mode.

Z	ion	I_i	compound	source	references
5	$^{11}\text{B}^{1+}$	4 nA*	$(\text{H}_3\text{CO})_3\text{B}$	EBIS	[Sasse 2012]
	$^{11}\text{B}^{1+}$	25 μA	$\text{C}_2\text{H}_{12}\text{B}_{10}$	ECRIS	[Bogomolov et al. 2015]
	$^{11}\text{B}^{2+}$	2 nA*	$(\text{H}_3\text{CO})_3\text{B}$	EBIS	[Sasse 2012]
	$^{11}\text{B}^{2+}$	70 μA	$\text{C}_2\text{H}_{12}\text{B}_{10}$	ECRIS	[Bogomolov et al. 2015]
	$^{11}\text{B}^{3+}$	1 nA*	$(\text{H}_3\text{CO})_3\text{B}$	EBIS	[Sasse 2012]
	$^{11}\text{B}^{3+}$	190 μA	$\text{C}_2\text{H}_{12}\text{B}_{10}$	ECRIS	[Bogomolov et al. 2015]
	$^{11}\text{B}^{4+}$	200 pA*	$(\text{H}_3\text{CO})_3\text{B}$	EBIS	[Sasse 2012]
	$^{11}\text{B}^{4+}$	85 μA	$\text{C}_2\text{H}_{12}\text{B}_{10}$	ECRIS	[Bogomolov et al. 2015]
	$^{11}\text{B}^{5+}$	60 pA*	$(\text{H}_3\text{CO})_3\text{B}$	EBIS	[Sasse 2012]
14	$^{28}\text{Si}^{5+}$	24 μA	$\text{Si}(\text{CH}_3)_4$	ECRIS	[Loginov et al. 2020]
	$^{28}\text{Si}^{5+}$	62 μA	$\text{Si}(\text{CH}_3)_4$	ECRIS	[Bondarchenko et al. 2020]
	$^{28}\text{Si}^{6+}$	40 μA	$\text{Si}(\text{CH}_3)_4$	ECRIS	[Bondarchenko et al. 2020]
	Si^{8+}	9 pA*	$\text{Si}[\text{Si}(\text{CH}_3)_3]_4$	EBIS	[Ullmann 2005]
	Si^{10+}	6 pA*	$\text{Si}[\text{Si}(\text{CH}_3)_3]_4$	EBIS	[Ullmann 2005]
	Si^{11+}	2 pA*	$\text{Si}[\text{Si}(\text{CH}_3)_3]_4$	EBIS	[Ullmann 2005]
	Si^{12+}	320 fA*	$\text{Si}[\text{Si}(\text{CH}_3)_3]_4$	EBIS	[Ullmann 2005]
	Si^{13+}	220 fA*	$\text{Si}[\text{Si}(\text{CH}_3)_3]_4$	EBIS	[Ullmann 2005]
	Si^{14+}	100 fA*	$\text{Si}[\text{Si}(\text{CH}_3)_3]_4$	EBIS	[Ullmann 2005]
15	P^{11+}	3 pA*	$(\text{C}_2\text{H}_5\text{O})_2\text{P}(\text{O})\text{H}$	EBIS	[Kentsch et al. 2004a]
	P^{15+}	255 fA*	$(\text{C}_2\text{H}_5\text{O})_2\text{P}(\text{O})\text{H}$	EBIS	[Kentsch et al. 2004a]

Z	ion	I_i	compound	source	references
22	$^{48}\text{Ti}^{5+}$	79 μA	$\text{C}_5(\text{CH}_3)_5\text{Ti}(\text{CH}_3)_3$	ECRIS	[Bogomolov et al 2014]
	$^{50}\text{Ti}^{5+}$	82 μA	$\text{C}_5(\text{CH}_3)_5\text{Ti}(\text{CH}_3)_3$	ECRIS	[Bogomolov et al 2014]
	$^{48}\text{Ti}^{9+}$	18 μA	$\text{C}_5(\text{CH}_3)_5\text{Ti}(\text{CH}_3)_3$	ECRIS	[Bondarchenko et al. 2020]
	$^{48}\text{Ti}^{10+}$	20 μA	$\text{C}_5(\text{CH}_3)_5\text{Ti}(\text{CH}_3)_3$	ECRIS	[Jardin et al. 2012]
	$^{48}\text{Ti}^{10+}$	16 μA	$\text{C}_5(\text{CH}_3)_5\text{Ti}(\text{CH}_3)_3$	ECRIS	[Bondarchenko et al. 2020]
	$^{48}\text{Ti}^{11+}$	68 μA	$\text{C}_5(\text{CH}_3)_5\text{Ti}(\text{CH}_3)_3$	ECRIS	[Bogomolov et al 2014]
23	V^{5+}	75 μA	$\text{V}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2016]
	V^{6+}	54 μA	$\text{V}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2016]
	V^{7+}	41 μA	$\text{V}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2016]
	V^{8+}	54 μA	$\text{V}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2016]
	V^{9+}	55 μA	$\text{V}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2016]
	V^{10+}	43 μA	$\text{V}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2016]
	V^{11+}	34 μA	$\text{V}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2016]
	V^{12+}	19 μA	$\text{V}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2016]
24	Cr^{5+}	50 μA	$\text{Cr}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2016]
	Cr^{6+}	70 μA	$\text{Cr}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2016]
	Cr^{7+}	60 μA	$\text{Cr}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2016]
	Cr^{8+}	37 μA	$\text{Cr}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2016]
	$^{52}\text{Cr}^{8+}$	22 μA	$\text{Cr}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2020]
	Cr^{9+}	17 μA	$\text{Cr}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2016]
	Cr^{10+}	7 μA	$\text{Cr}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2016]
	$^{52}\text{Cr}^{10+}$	9 μA	$\text{Cr}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2020]
	$^{52}\text{Cr}^{10+}$	18 μA	$\text{Cr}(\text{C}_5\text{H}_5)_2$	ECRIS	[Loginov et al. 2020]
26	Fe^{6+}	43 μA	$\text{Fe}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2016]
	Fe^{6+}	7 μA^{lm}	$\text{Fe}(\text{C}_5\text{H}_5)_2$	EBIS	[Sasse 2012]
	Fe^{7+}	93 μA	$\text{Fe}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2016]
	Fe^{8+}	125 μA	$\text{Fe}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2016]
	$^{56}\text{Fe}^{8+}$	68 μA	$\text{Fe}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2020]
	Fe^{9+}	172 μA	$\text{Fe}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2016]
	Fe^{10+}	145 μA	$\text{Fe}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2016]
	$^{56}\text{Fe}^{10+}$	44 μA	$\text{Fe}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2020]

Z	ion	I_i	compound	source	references
	Fe ¹⁰⁺	2 μA^{lm}	Fe(C ₅ H ₅) ₂	EBIS	[Sasse 2012]
	Fe ¹¹⁺	114 μA	Fe(C ₅ H ₅) ₂	ECRIS	[Bondarchenko et al. 2016]
	⁵⁶ Fe ¹¹⁺	25 μA	Fe(C ₅ H ₅) ₂	ECRIS	[Bondarchenko et al. 2020]
	Fe ¹²⁺	73 μA	Fe(C ₅ H ₅) ₂	ECRIS	[Bondarchenko et al. 2016]
	⁵⁶ Fe ¹²⁺	12 μA	Fe(C ₅ H ₅) ₂	ECRIS	[Bondarchenko et al. 2020]
	Fe ¹²⁺	800 pA ^{lm}	Fe(C ₅ H ₅) ₂	EBIS	[Sasse 2012]
	Fe ¹³⁺	45 μA	Fe(C ₅ H ₅) ₂	ECRIS	[Bondarchenko et al. 2016]
	Fe ¹⁷⁺	20 pA*	Fe(C ₅ H ₅) ₂	EBIS	[Sasse 2012]
	Fe ¹⁷⁺	300 pA ^{lm}	Fe(C ₅ H ₅) ₂	EBIS	[Sasse 2012]
	Fe ²⁰⁺	100 pA ^{lm}	Fe(C ₅ H ₅) ₂	EBIS	[Sasse 2012]
	Fe ²⁴⁺	5 pA ^{lm}	Fe(C ₅ H ₅) ₂	EBIS	[Sasse 2012]
	Fe ²⁵⁺	200 fA ^{lm}	Fe(C ₅ H ₅) ₂	EBIS	[Sasse 2012]
	Fe ²⁵⁺	625 fA*	Fe(C ₅ H ₅) ₂	EBIS	[Sasse 2012]
	Fe ²⁶⁺	520 fA*	Fe(C ₅ H ₅) ₂	EBIS	[Sasse 2012]
27	Co ⁶⁺	57 μA	Co(C ₅ H ₅) ₂	ECRIS	[Bondarchenko et al. 2016]
	Co ⁷⁺	80 μA	Co(C ₅ H ₅) ₂	ECRIS	[Bondarchenko et al. 2016]
	Co ⁸⁺	86 μA	Co(C ₅ H ₅) ₂	ECRIS	[Bondarchenko et al. 2016]
	Co ⁹⁺	98 μA	Co(C ₅ H ₅) ₂	ECRIS	[Bondarchenko et al. 2016]
	⁵⁹ Co ⁹⁺	57 μA	Co(C ₅ H ₅) ₂	ECRIS	[Bondarchenko et al. 2020]
	⁵⁹ Co ¹⁰⁺	73 μA	Co(C ₅ H ₅) ₂	ECRIS	[Bondarchenko et al. 2020]
	Co ¹¹⁺	82 μA	Co(C ₅ H ₅) ₂	ECRIS	[Bondarchenko et al. 2016]
	⁵⁹ Co ¹¹⁺	36 μA	Co(C ₅ H ₅) ₂	ECRIS	[Bondarchenko et al. 2020]
	Co ¹²⁺	25 μA	Co(C ₅ H ₅) ₂	ECRIS	[Bondarchenko et al. 2016]
	⁵⁹ Co ¹²⁺	12 μA	Co(C ₅ H ₅) ₂	ECRIS	[Bondarchenko et al. 2020]
	⁵⁹ Co ¹²⁺	12 μA	Co(C ₅ H ₅) ₂	ECRIS	[Loginov et al. 2020]
28	Ni ⁶⁺	45 μA	Ni(C ₅ H ₅) ₂	ECRIS	[Bondarchenko et al. 2016]
	Ni ⁷⁺	43 μA	Ni(C ₅ H ₅) ₂	ECRIS	[Bondarchenko et al. 2016]
	Ni ⁸⁺	48 μA	Ni(C ₅ H ₅) ₂	ECRIS	[Bondarchenko et al. 2016]
	Ni ⁹⁺	53 μA	Ni(C ₅ H ₅) ₂	ECRIS	[Bondarchenko et al. 2016]
	⁵⁸ Ni ⁹⁺	39 μA	Ni(C ₅ H ₅) ₂	ECRIS	[Bondarchenko et al. 2020]
	Ni ¹¹⁺	30 μA	Ni(C ₅ H ₅) ₂	ECRIS	[Bondarchenko et al. 2016]

Z	ion	I_i	compound	source	references
	$^{58}\text{Ni}^{11+}$	43 μA	$\text{Ni}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2020]
	$^{58}\text{Ni}^{11+}$	43 μA	$\text{Ni}(\text{C}_5\text{H}_5)_2$	ECRIS	[Loginov et al. 2020]
	Ni^{12+}	10 μA	$\text{Ni}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2016]
	$^{58}\text{Ni}^{12+}$	37 μA	$\text{Ni}(\text{C}_5\text{H}_5)_2$	ECRIS	[Bondarchenko et al. 2020]
	$^{58}\text{Ni}^{12+}$	9 μA	$\text{Ni}(\text{C}_5\text{H}_5)_2$	ECRIS	[Koivisto et al. 1998]
	$^{58}\text{Ni}^{13+}$	6 μA	$\text{Ni}(\text{C}_5\text{H}_5)_2$	ECRIS	[Koivisto et al. 1998]
	$^{58}\text{Ni}^{14+}$	3 μA	$\text{Ni}(\text{C}_5\text{H}_5)_2$	ECRIS	[Koivisto et al. 1998]
	$^{58}\text{Ni}^{15+}$	1 μA	$\text{Ni}(\text{C}_5\text{H}_5)_2$	ECRIS	[Koivisto et al. 1998]
32	Ge^{7+}	43 μA	$\text{Ge}(\text{CH}_2\text{CH}_3)_4$ $\text{Ge}(\text{CH}_3)_4$	ECRIS	[Bondarchenko et al. 2016]
	Ge^{8+}	54 μA	$\text{Ge}(\text{CH}_2\text{CH}_3)_4$ $\text{Ge}(\text{CH}_3)_4$	ECRIS	[Bondarchenko et al. 2016]
	$^{72}\text{Ge}^{8+}$	38 μA	$\text{Ge}(\text{CH}_3)_4$	ECRIS	[Bondarchenko et al. 2020]
	Ge^{10+}	47 μA	$\text{Ge}(\text{CH}_2\text{CH}_3)_4$ $\text{Ge}(\text{CH}_3)_4$	ECRIS	[Bondarchenko et al. 2016]
	$^{72}\text{Ge}^{10+}$	46 μA	$\text{Ge}(\text{CH}_3)_4$	ECRIS	[Bondarchenko et al. 2020]
	$^{72}\text{Ge}^{11+}$	34 μA	$\text{Ge}(\text{CH}_3)_4$	ECRIS	[Bondarchenko et al. 2020]
72	Hf^{13+}	31 μA	$(\text{C}_5\text{H}_5)_2\text{Hf}(\text{CH}_3)_2$	ECRIS	[Bondarchenko et al. 2016]
	Hf^{14+}	45 μA	$(\text{C}_5\text{H}_5)_2\text{Hf}(\text{CH}_3)_2$	ECRIS	[Bondarchenko et al. 2016]
	Hf^{16+}	50 μA	$(\text{C}_5\text{H}_5)_2\text{Hf}(\text{CH}_3)_2$	ECRIS	[Bondarchenko et al. 2016]
	Hf^{17+}	45 μA	$(\text{C}_5\text{H}_5)_2\text{Hf}(\text{CH}_3)_2$	ECRIS	[Bondarchenko et al. 2016]
	$^{180}\text{Hf}^{17+}$	19 μA	$(\text{C}_5\text{H}_5)_2\text{Hf}(\text{CH}_3)_2$	ECRIS	[Bondarchenko et al. 2020]
	Hf^{18+}	36 μA	$(\text{C}_5\text{H}_5)_2\text{Hf}(\text{CH}_3)_2$	ECRIS	[Bondarchenko et al. 2016]
	$^{180}\text{Hf}^{18+}$	15 μA	$(\text{C}_5\text{H}_5)_2\text{Hf}(\text{CH}_3)_2$	ECRIS	[Bondarchenko et al. 2020]
	Hf^{19+}	27 μA	$(\text{C}_5\text{H}_5)_2\text{Hf}(\text{CH}_3)_2$	ECRIS	[Bondarchenko et al. 2016]
	Hf^{20+}	17 μA	$(\text{C}_5\text{H}_5)_2\text{Hf}(\text{CH}_3)_2$	ECRIS	[Bondarchenko et al. 2016]
	$^{180}\text{Hf}^{20+}$	5 μA	$(\text{C}_5\text{H}_5)_2\text{Hf}(\text{CH}_3)_2$	ECRIS	[Bondarchenko et al. 2020]
80	$^{202}\text{Hg}^{24+}$	20 pA*	Hg	EBIS	[Sasse 2012]
	$^{202}\text{Hg}^{33+}$	9 pA*	Hg	EBIS	[Sasse 2012]
	$^{202}\text{Hg}^{43+}$	4 pA*	Hg	EBIS	[Sasse 2012]
	$^{202}\text{Hg}^{51+}$	1 pA*	Hg	EBIS	[Sasse 2012]
	$^{202}\text{Hg}^{55+}$	152 fA*	Hg	EBIS	[Sasse 2012]

Z	ion	I_i	compound	source	references
	Hg ⁵⁸⁺	70 fA*	Hg	EBIS	[Sasse 2012]

References

- [Ärje et al. 1993] J. Ärje, H. Koivisto, M. Nurmi, Proceedings of the 11th Workshop on ECR Ion Sources, KVI-Report 996 (1993) 27
- [Barue et al. 2000] C. Barue, C. Canet, M. Dupuis, J.L. Flambard, P.A. Jaffres, P. Leherissier, F. Lemagnen, Workshop on the Production of Intense Beams of Highly Charged Ions, September 24-27, 2000, Catania, Italy, GANIL P 00 42
- [Bogomolov et al. 1999] S.L. Bogomolov, A.A. Efremov, V.B. Kutner, A.N. Lebedev, V.N. Loginov, N. Yu, *Production of ions of metals with an ECR ion source at FLNR (JINR) cyclotrons*, in: Proceedings of the 14th International Workshop on ECR Ion Sources, CERN, Switzerland, May 3-6, 1999, 71
- [Bogomolov et al. 2014] S. Bogomolov, A. Bondarchenko, A. Efremov, K. Kuzmenkov, A. Lebedev, K. Lebedev, V. Lebedev, V. Loginov, N. Yazvitsky, Z. Asfari, B.J.P. Gall, *Production of metal ion beams from ECR ion sources by MIVOC method*, XXIV Russian Particle Accelerator Conference RuPAC 2014, 6-10 October 2014, Obninsk, Russia, 432 (THPSC47)
- [Bogomolov et al. 2015] S.L. Bogomolov, A.E. Bondarchenko, A.A. Efremov, K.I. Kuzmenkov, A.N. Lebedev, K.V. Lebedev, V.Ya. Lebedev, V.N. Loginov, V.E. Mironov, N.Yu. Yazvitsky, *Physics of Particles and Nuclei Letters* 12, No. 7 (2015) 824
- [Bondarchenko et al. 2016] A.E. Bondarchenko, S.L. Bogomolov, K.I. Kuzmenkov, V.N. Loginov, V.Ya. Lebedev, Z. Asfari, B.J.P. Gall, Proceedings of the Russian Particle Accelerator Conference (RuPAC'16), St. Petersburg, Russia, November 21-25, 2016, TUPSA059, 330
- [Bondarchenko et al. 2020] A. Bondarchenko, S. Bogomolov, V. Loginov, A. Lebedev, V. Mironov, D. Pugachev, M. Zdorovets, I. Ivanov, Y. Sambyev, M. Koloberdin, A. Kurakhmedov, D. Mustafin, M. Abdigaliyev, Proceedings of the 24th International Scientific Conference of Young Scientists and Specialists (AYSS-2020), AIP Conf. Proc. 2377, 080001-4, <https://doi.org/10.1063/5.0063412>
- [Jardin et al. 2012] P. Jardin et al., *"ECRISs at GANIL today and tomorrow"*, in Proceedings of ECRIS2012, Sydney, Australia, p.195; <http://accelconf.web.cern.ch/Accel-Conf/ECRIS2012/papers/frya01.pdf>
- [Jovovic et al. 2007] J. Jovovic, J. Cvetic, A. Dobrosavljevic, T. Nedeljkovic, I. Draganic, Proc. 51st ETRAN Conference, Herceg Novi-Igalo, June 4-8, 2007
- [Kentsch et al. 2002] U. Kentsch, G. Zschornack, F. Grossmann, V.P. Ovsyannikov, F. Ullmann, *Review of Scientific Instruments* 73 (2002) 657
- [Kentsch et al. 2002a] U. Kentsch, G. Zschornack, F. Grossmann, V.P. Ovsyannikov, F. Ullmann, S. Fritzsche, A. Surzhykov, *Nuclear Instruments & Methods B187* (2002) 238

- [Kentsch et al. 2004] U. Kentsch, S. Landgraf, G.Zschornack, F. Grossmann, V.P. Ovsyannikov, F. Ullmann, X-Ray Spectrometry 33 (2004) 33
- [Kentsch et al. 2004a] U.Kentsch, S. Landgraf, M. Schmidt, H. Tyrroff, G. Zschornack, F. Grossmann, V.P. Ovsyannikov, F. Ullmann, Nuclear Instruments & Methods B216 (2004) 196
- [Koivisto et al. 1994] H. Koivisto, J. Ärje, M. Nurmia, Nuclear instruments & Methods in Physics research B94 (1994) 291
- [Koivisto et al. 1998] H.Koivisto, J.Ärje, M. Nurmia, Review of Scientific Instruments 69 (1998) 785
- [Koivisto et al. 2001] H. Koivisto, P. Heikkinen, V. Hänninen, A. Lassila, H. Leinonen, V. Nieminen, J. Pakarinen, K. Ranttila, J. Ärje, E. Liukkonen, Nuclear Instruments & Methods in Physics Research B174 (2001) 379
- [Koivisto et al. 2002] H. Koivisto, J. Arje, R. Seppala, M. Nurmia, Nuclear Instruments & Methods in Physics Research B187 (2002) 111
- [Liu et al 2020] H.P.Liu, J.Y. Li, L. Kang, T.M. Zhang, J. Li, IOP Conference Series: Materials Science and Engineering 733 (2020) 012034, doi:1088/1757-899X/733/1/012034
- [Loginov et al. 2018] V.N. Loginov, A.E. Bondarchenko, K.I. Kuzmenkov, S.L. Bogomolov, V.E. Mironov, I.A. Ivanov, M.V. Zdorovets, E.K. Sambayev, V.V. Alexandrenko, A.E. Kurakhmedov, S.G. Kozin, M.V. Koloberdin, D.A. Mustafin, Proceedings of the 26th Particle Accelerator Conference RUPAC 2018, Protvino, Russia, p. 504 (THPSC50), doi:10.18429/JACoW-RUPAC2018-THPSC50
- [Loginov et al. 2020] V.N. Loginov, S.L. Bogomolov, A.E. Bondarchenko, V.E. Mironov, D.K. Pugachev, The Journal Physics of Elementary Particles and Atomic Nuclei, Letters, No. 17, Issue 2 (2020) 153
- [Lyman et al. 1990] W.J. Lyman, W.F. Reehl, D.H. Rosenblatt, Handbook of Chemical Property Estimation Methods, American Chemical Society, Washington, DC, McGraw-Hill 1990
- [Nakagawa et al. 1998] T. Nakagawa, J. Ärje, Y. Miyazawa, M. Hemmi, T. Chiba, N. Inabe, M. Kase, T. Kageyama, O. Kamigaito, M. Kidera, A. Got, Y. Yano, Review of Scientific Instruments 69 (1998) 637
- [Nakamura et al. 2000] N. Nakamura, T. Kinugawa, H. Shimizu, H. Watanabe, S. Ito, S. Ohtani, C. Yamada, K. Okazaki, M. Sakurai, M.R. Tarbutt, J.D. Silver, Review of Scientific Instruments 71 (2000) 684
- [Oshiro et al. 2014] Y. Oshiro, S. Yamaka, S. Watanabe, K. Kobayashi, Y. Kotaka, M. Nishimura, M. Kase, H. Muto, H. Yamaguchi, S. Shimoura, Review of Scientific Instruments 85 (2014) 02A912

- [Pugachev et al. 2021] D.K. Pugachev, S.L. Bogomolov, V.E. Mironov, A.A. Efremov, V.N. Loginov, A.N. Lebedev, A.E. Bondarchenko, K.I. Kuzmmenkov, K.B. Gikal, A.A. Protasov, 27th Russian Particle Accelerator Conference RuPAC2021, Alushta, Russia, 2021, TUPSB36, p. 303, doi:10.18429/JACoW-RuPAC2021-TUPSB36
- [Sasse 2012] H.Sasse, *Untersuchungen zu Beladungstechniken von EBIS mit Metallionen*, TU Dresden, Fakultät Mathematik und Naturwissenschaften, Fachrichtung Physik, Dresden, 2012
- [Stiebing et al. 1999] K.E. Stiebing, O.Hohn, S.Runkel, L. Schmidt, H. Schmidt-Böcking, V. Mironov, G. Shirkov, Physical Review Special Topics - Accelerators and Beams 2 (1999) 123501
- [Takasugi et al. 2010] W. Takasugi, M. Wakaisami, N. Sasaki, T. Sakuma, M. Yamamoto, A. Kitagawa, M. Muramatsu, Review of Scientific instruments 81 (2010) 02A329
- [Ullmann 2005] F. Ullmann, *Untersuchung der Erzeugung hochgeladener Ionen in einer Raumtemperatur-Elektronenstrahl-Ionenfalle*, Dissertation, TU Dresden, 2005
- [Ullrich 1997] A. Ullrich, *Wellenlängendispersive Röntgenspektroskopie am mit Metallionen beladenen ECR-Plasma zur Bestimmung der Ionenladungsverteilung im Plasma*, Diploma Thesis, TU Dresden, Fachrichtung Physik, 1997
- [Ullrich et al. 1998] A.Ullrich, P. Grübling, G. Zschornack, Review of Scientific Instruments 69 (1998) 813
- [Weast 1978] R.C. Weast, Handbook of Chemistry and Physics, 58th Edition, CRC Press, 1978 p. D-176
- [Werner et al. 2000] T. Werner, G. Zschornack, G. Grossmann, V.P. Ovsyannikov, F.Ullmann, Nuclear Instruments & Methods B178 (2000) 2038
- [ZiMin et al. 2000] Z. ZiMin, Z. HongWei, Z. YuBin, Z. XueZhen, G. Xiaohong, L. ZhanWen, Chinese Physics C 24 (2000) 974